

Lab#6

Microfabrication Principles

Micromolding and Softlithography

1. Check Lab Notebooks and Collect Project Assignments.
2. This lab has two parts, which are 1. Fabrication of mold using SU-8 and 2. Fabrication of microstructures using Softlithography.
3. You will use one polished silicon wafer per batch.
4. Part#1 Fabrication of SU-8 Mold
 - i. Remove SU-8 10 from the refrigerator for at least 2 hours before starting the process.
 - ii. Start with cleaning process. (Refer lab#1/lab#5 cleaning procedure)
 - iii. Bake the silicon wafer at 200°C for 5 minutes.
 - iv. Dispense 3 ml SU-8 10 on the wafer (3 inch) after placing the silicon wafer on the spinner.
 - v. Spin at (I) 200 rpm, 10 seconds and stop. (II) 1000 rpm, 30 seconds. This spin speed will provide ~ 25 μm thickness.
 - vi. Soft bake at 65°C for 3 minutes and 90°C for 7 minutes.
 - vii. Place the wafer on 50°C hot plate for 2 minutes and then allow it to cool to room temperature.
 - viii. Exposure dose is 230 -250 mJ/cm^2 . Use the optical power meter to measure the lamp intensity and figure out the exposure time.
 - ix. Post Bake (I) 65°C for 1 minute and (II) 95°C for 3 minute.
 - x. Develop in SU-8 developer for 4 minutes.
 - xi. Clean the SU-8 developer with IPA (avoid acetone.), and then dry with air or nitrogen.
 - xii. Hard bake at 150°C for 15 minutes.
5. Part#2 Fabrication of PDMS microstructures using soft lithography.
 - i. Add PDMS monomer and curing agent in 10:1 by weight proportion and stir vigorously. We will do 10g base and 1g curing agent.
 - ii. Keep the mixture in the vacuum chamber for 30 minutes to remove air bubbles.
 - iii. Pour degassed elastomer mixture on a clean and polished silicon wafer.
 - iv. Carefully place the wafer in oven at 150°C for 10 min or until PDMS cures.
6. Slowly peel off PDMS microstructure. Measure and compare the height of the PDMS and SU-8 microstructures.